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LABEL TEXT

PREDNISONE TABLETS 30 tablets

The nominal weight of prednisone in each tablet is 10 mg. Use only whole tablets. See the certificate for additional information. Store in a dry place. Store at controlled room temperature not exceeding 25°C. Unused or unopened blister strips should be kept in the secondary package.

Danger! Causes eye irritation. Suspected of damaging fertility or the unborn child. Causes damage to organs (endocrine system) through prolonged or repeated exposure.

Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Wash thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. If exposed or concerned: Get medical advice/attention. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations.

For use with specified USP compendial tests. Not for use as a drug. See SDS prior to use at www.usp.org/ods.

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Apparatus	# of vessels	Single-Stage		Two-Stage			
				1 st Stage of Two Stages		2 nd Stage of Two Stages	
		GM*	%CV	GM*	%CV	GM*	%CV
1	6	56-88	12	61-81	8.6	56-88	11
	7		11			56-87	
	8	56-87					
	12	na					
	14	56-88		na			
2	6	28-39	7.3	29-37	5.5	28-39	7.2
	7		7.2			7.1	
	8		7.0			7.0	
	12		7.2	na			
	14		7.1	na			

* Percent of the labeled amount of prednisone dissolved at 30 minutes at 50-rpm

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1st Stage of Two-Stage for $n=6, 7, 8$ and Single-Stage for $n=12, 14$:
 $GM1 = \exp(\text{average}(\ln x_1 : \ln x_n))$
 $\%CV1 = 100 * \sqrt{\exp(\text{var}(\ln x_1 : \ln x_n)) - 1}$

Single-Stage or 2nd Stage of Two-Stage for $n=6, 7, 8$:
 $GM = \exp(\text{average}((\text{average}(\ln x_1 : \ln x_n)), (\text{average}(\ln x_{n+1} : \ln x_{2n})))) = \exp(\text{average}(\ln x_1 : \ln x_{2n}))$
 $\%CV = 100 * \sqrt{\exp(\text{average}((\text{var}(\ln x_1 : \ln x_n)), (\text{var}(\ln x_{n+1} : \ln x_{2n})))) - 1}$

exp: exponential (e^x) var: variance sqrt: square root *: multiply 100: conversion factor to percentage

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Performance Verification Test Calculation Tool

Q1. Please select the Reference Standard Lot: R05991 ▼

Q2. Which Dissolution Apparatus are you using?

Basket
Paddle

Q3. Which configuration of dissolution test assemblies are you using?

6
7
8
12
14

Q4. Please choose the testing procedure you will be using to perform your calculation?

1-stage
2-stage

Submit Data

PVT Calculation Tool V2.0

Apparatus 1 -- 1-stage test for equipment with 6 positions

This worksheet is for use for the PVT for Apparatus 1 with 6 positions where a decision has been made to use the 1-stage option

The valid use of the PVT requires that the choice of 1- or 2-stage testing is made before beginning testing. The use of this worksheet is not valid if the choice of 1- or 2-stage testing is made after examining the data.

Acceptance criteria for USP Prednisone Tablets RS, R05991

GeoMean, lower	56
GeoMean, upper	88
%CV, upper	12

Enter data as amount dissolved at 30 minutes as percent of 10 mg.
(Data entered with more than three decimals will be rounded to three.)

Run 1 Data

Vessel 1

Vessel 2

Vessel 3

Vessel 4

Vessel 5

Vessel 6

Run 2 Data

Vessel 1

Vessel 2

Vessel 3

Vessel 4

Vessel 5

Vessel 6

Submit Data

Apparatus 2 -- 1-stage test for equipment with 6 positions

This worksheet is for use for the PVT for Apparatus 2 with 6 positions where a decision has been made to use the 1-stage option

The valid use of the PVT requires that the choice of 1- or 2-stage testing is made before beginning testing. The use of this worksheet is not valid if the choice of 1- or 2-stage testing is made after examining the data.

Acceptance criteria for USP Prednisone Tablets RS, R05991

GeoMean, lower	28
GeoMean, upper	39
%CV, upper	7.3

Enter data as amount dissolved at 30 minutes as percent of 10 mg.
(Data entered with more than three decimals will be rounded to three.)

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1. United States Pharmacopeial Convention. <711> Dissolution. The United States Pharmacopeia - National Formulary: USP 40 - NF 35 2017. 35th ed. Baltimore, MD: The United Book Press, Inc.; 2016. p. 588-98.
2. United States Pharmacopeial Convention. <1058> Analytical Instrument Qualification. The United States Pharmacopeia - National Formulary: USP 40 - NF 35 2017. 35th ed. Baltimore, MD: The United Book Press, Inc.; 2016. p. 1137-43.
3. Dissolution Toolkit Procedures for Mechanical Calibration and Performance Verification Test Apparatus 1 and Apparatus 2. [Internet]. 2010 [cited 2017 Dec 5]. Available from: <http://www.usp.org/sites/default/files/usp/document/our-work/reference-standards/dissolution-toolkit-version2.pdf>.
4. Certificate of Prednisone Tablets. [Internet]. 2017 [cited 2017 Dec 12]. Available from: http://store.usp.org/OA_HTML/usp3_ibeuspcertiPrevLot.jsp?&CatalogNo=1559505.
5. Calculation Tool for Performance Verification Test Calculations. [Internet]. 2017 [cited 2017 Dec 12]. Available from: <http://apps.usp.org/app/USPNF/pvtCalculationTool/>.